

Implementing Study Amendments Without Downtime or Data Migrations

Transforming database design with a modern EDC

The Challenge of Migrations in Clinical Trials

The growing complexity in clinical trials is driving an increase in protocol amendments. Research from the Tufts Center for the Study of Drug Development found that Phase III studies have an average of 3.6 protocol amendments, each containing 8.5 changes.¹ That is over 30 individual changes in total. This count does not include other changes, such as corrections, updates to procedures, and new edit checks, meaning the total number of structural and functional changes is much higher.

For most studies, those amendments and other changes all required migrations.

An article in Applied Clinical Trials describing [how study complexity is directly correlated with the number of changes needed](#) cited that studies with simple Case Report Forms (CRFs) required an average of eight migrations and those with highly complex CRFs saw an average of 17 migrations.² These numbers are concerning when you consider the impact of changes and migrations on study timelines. According to one industry whitepaper, “The typical turn-around time for implementing a limited change, for example, adding 1 existing form in a new visit, is 10 working days.”³ Certainly the time required to implement the 8.5 changes in a typical amendment is much higher.

¹ Rising Protocol Complexity is Hindering Performance while Driving Up Cost of Clinical Trials, Tufts CSDD, Impact Report July/August 2018, Vol. 20 No.4.

² Time for Another CRF Migration, Medidata Solutions, Applied Clinical Trials, December 2014. Accessed via: <http://www.appliedclinicaltrialsonline.com/time-another-crf-migration>

³ “Innovative approach to building an adaptive trial design in Medidata Rave®” Joris De Bondt, SGS Life Science Services, Paper TS07, PhUSE 2015.

Why Changes Required Time-consuming Migrations

Implementing changes in traditional EDC systems takes a long time because organizations must build a new database and physically move existing clinical data. To achieve an effective data migration, data fields are mapped between systems, data is extracted from the old system and loaded to the new system.

After loading into the new system, results must be verified to determine whether data was accurately translated, is complete, and supports processes in the new system. During verification, you may need to run both systems in parallel to prevent accidental data loss while checking for disparities.

The system must be down or users blocked for the duration of all migration and test procedures. Thus, each migration introduces study downtime causing delays that can sometimes extend over several days.

The authors in the above mentioned ACT article [advised companies to simplify their eCRFs in order to minimize the burden of migration costs and downtime](#). However, dumbing-down eCRFs creates risks around data quality, sub-optimal study designs, and poor alignment to increasingly complex protocols. This situation underlies the need for a better EDC, one that not only simplifies change management but embraces it.

A Modern Architecture and Approach in Veeva Vault EDC

A modern EDC allows for a casebook design that is as dynamic as the clinical trial itself. Veeva has delivered a next-generation EDC system, designed specifically to address the migration issue. Data is stored in a manner that provides flexibility and removes the need for migrations entirely. The ability to insert new data collection elements, new pages, and new visits, without bringing the system down or migrating data, simplifies the initial build as well as subsequent amendments.

A Flexible Architecture that Supports Change

With Veeva, all visits, forms, and fields are defined using a single data model with database tables and columns that do not change between studies. New studies and their eCRFs are defined by adding rows to the existing tables and columns. Once the study is live, you can add or change elements by adding new rows, without changing the database structure.


Vault EDC uses reference IDs as keys to the underlying data rather than using the physical location of the data. When clinical data is entered, it is stored with a reference pointer to the definition for that field. When a field is changed, a new reference is created for the stored value. As a result, the system can have multiple versions of a casebook at any time, there are no changes to the database structure, and study migrations are obviated.

Deploying Changes Without Disruption

Vault EDC saves you time and money by implementing changes in near real-time. When deploying a change, Vault automatically locks users out of a casebook for the few minutes it takes to effect the change and then releases the page. Normal service is resumed within minutes. And by scheduling deployments during non-working hours at each site, there is zero disruption for sites.

The benefits from this flexibility are considerable: the desired changes go into effect faster, administration costs are minimal, and downtime for sites is eliminated. Perhaps most importantly, your trials become more agile and making amendments no longer poses risk to your data.

Figure 3. Multiple study design versions coexist within the database and are assigned with patient-level granularity.



About Veeva Systems

Veeva is the global leader in cloud software for the life sciences industry. Committed to innovation, product excellence, and customer success, Veeva serves more than 1,100 customers, ranging from the world's largest pharmaceutical companies to emerging biotechs. As a Public Benefit Corporation, Veeva is committed to balancing the interests of all stakeholders, including customers, employees, shareholders, and the industries it serves. For more information, visit www.veeva.com.

Veeva Systems

Global Headquarters
Pleasanton, California, USA
4280 Hacienda Drive
Pleasanton, California 94588
+1 925 452 6500 | veeva.com/contact | veeva.com

Copyright © 2026 Veeva Systems. All rights reserved. Veeva and the Veeva logo are registered trademarks of Veeva Systems. Veeva Systems owns other registered and unregistered trademarks. Other names used herein may be trademarks of their respective owners.

Learn more at veeva.com | 925-452-6500 | veeva.com/contact